

# Mobile X-Ray System for Infection Control



MobileDaRt Evolution  
MX8<sup>Version</sup>

The MobileDaRt MX8 is your optimal solution for enabling X-Ray imaging in any emergency

## Benefits of the MobileDaRt Evolution MX8 with Digital Wireless Panel

- Display images on the large 19" monitor within 2 seconds after exposure
- Images can be sent to a physician at a remote location for quick diagnosis (both wired and wireless connection require hospital internet connection)
- System boot-up takes only 1 minute and is immediately ready for use in emergency conditions
- Easily and quickly disinfected with common hospital alcohols and cleaners\*

\*Reference Shimadzu document M516-E056 (Cleaning Instructions to Avoid Secondary Infection)



## Contact Shimadzu Medical Systems USA

20101 S. Vermont Avenue - Torrance, CA Phone: 800-228-1429  
[COVID19SPECIALRESPONSE@shimadzu-usa.com](mailto:COVID19SPECIALRESPONSE@shimadzu-usa.com)






Reliability and confidence from  
Shimadzu Medical Systems

## MobileDaRt Evolution MX8 Radiography System



During a pandemic infection, it is expected that the existing X-ray room will not be available to maintain isolation between infected and uninfected people, therefore it is necessary to prepare new isolated area(s) inside the hospital or outside in tented controlled spaces, and to have an X-Ray system that can be operated in that area.



### Easy Charge Management

The MobileDaRt MX8 is powered by built-in batteries (FPD needs to be charged separately from the main body). In an isolated area, you can easily accomplish the complete workflow exposure to image confirmation and data transfer all in one unit (requires hospital internet connection) The Mobile X-ray system can continue to be used even while charging.

email:

[COVID19SPECIALRESPONSE@shimadzu-usa.com](mailto:COVID19SPECIALRESPONSE@shimadzu-usa.com)

