



Excellence in Science

# Shimadzu Model X Promotion



FLUOROsPEED X1 edition

SONIALVISION G4 LX edition



Shimadzu Medical Systems USA newly announces Shimadzu XLT panel release. Please experience this versatile, durable and reliable DR system and reduce your cost of ownership while improving productivity.

## XLT panel highlights

- 14"x17" versatile panel
- Light panel in its class 5.5lbs. (2.5kg)
- The flexibility and strength to hold up to 661lbs. on a bed
- Liquid resistance, IPX6 compliant
- 13min. of charge provides 4.1 hours of use

**SPECIAL OFFER**

**Limited offer until the end of March 2022**

**Total 2 years DR panel protection  
1 year license of system productivity analysis tool**

Contact us ([marketing@shimadzu-usa.com](mailto:marketing@shimadzu-usa.com)) to get more details.

Register to learn more about Shimadzu Healthcare Solutions  
<https://www.shimadzuhealthcaresolutions.com/dollhouse-view/>

Shimadzu Medical Systems USA

20101 S. Vermont Avenue

Torrance, CA 90502 (800) 228-1429

[www.shimadzu-usa.com](http://www.shimadzu-usa.com)

Email : [Email: marketing@shimadzu-usa.com](mailto:marketing@shimadzu-usa.com)



**We Care.  
For You.**

Technology innovation is at our core at Shimadzu, and we use it to solve problems in healthcare. For over 100 years of history in X-ray, we have pursued our passion for technology to develop solutions to lower dose, secure patient safety, make workflow simpler and improve the patient and healthcare worker experience with meaningful innovation. At Shimadzu, we will continue to challenge ourselves to apply technology for a more meaningful experience for healthcare system.



 **SHIMADZU**  
Excellence in Science

**Shimadzu Medical Systems USA**  
20101 S. Vermont Avenue  
Torrance, CA 90502 (800) 228-1429  
[www.shimadzu-usa.com](http://www.shimadzu-usa.com)  
Email : [marketing@shimadzu-usa.com](mailto:marketing@shimadzu-usa.com)